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29 SEP 2011

26 September 2011

**Ministry of Health and Social Development
Department of High-Technological Medical Aid
Deputy of Director
Mrs. Gabbasova L.A.**

Dear Mrs. Gabbasova

Cochlear Limited company (Australia), world leading manufacturer of cochlear implantation systems, informed you that it is undertaking a voluntary recall of any unimplanted stock of Cochlear Nucleus CI500 series cochlear implants.

At present 0,6% of Nucleus CI 500 series implants have failed. This percentage is the lowest among all rates of different implants by different manufacturers in the world market.

However Cochlear being serious and responsible manufacturer, helping more than 250 000 hearing impaired people in the world, being the first Company in cochlear implants research, never having any kind of recalls in its history, in an abundance of caution and with its recipients in mind organized a voluntary recall of CI500 series implants for further investigation of causes of this issue.

It is important to note that if a Nucleus CI500 series implant failure occurs, the implant safely shuts down. In the event of such a failure there is no negative effect for recipients' health.

CI500 series implants include: CI512 implants and also CI513, CI551, ABI541 which are not widely spread and used in a limited quantity of countries.

Implants outside of the Nucleus CI500 series are not affected by this voluntary recall: these are the Nucleus Freedom Series with Contour Advance™ CI24 RE(CA) or Straight Electrode CI24RE(ST), Nucleus Hybrid™ L24 CI24REH, the Nucleus CI422 Slim straight electrode, the Nucleus 24 R series CI24RCA, CI24RST, CI24RCS, the Nucleus 24M Double Array CI11+11+2M and the Nucleus ABI24M. No sound processors or the external devices (including CP810 sound processor of Nucleus 5 generation) are affected by this recall.

All existing recipients with a Nucleus CI500 Series implant can continue to use their system as normal.

At present Cochlear offers to the new recipients – candidates for the cochlear implantation – the CI system in the following optimal configuration: CI24RE series cochlear implant (internal part) of reliability, sound processor Nucleus 5 series - CP810 (external part).

CI24RE Series implant shows extremely positive statistics of reliability and outstanding hearing performance results. To date, more than 70'000 recipients are now using the CI24RE series implant.

For more than 25 years cochlear implant reliability has been of paramount importance to Cochlear and we remain committed to providing high quality, reliable cochlear implant systems. All necessary measures are now being undertaken for eliminating the causes of failures of CI500 Series implants and for further improvement in its reliability statistics.

Yours faithfully,

**Richard Brook
President Europe, Middle East and Africa
Cochlear AG, Cochlear Limited**

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